BIOLOGIE CELULARA SI MOLECULARA

LP06. Izolare acizi nucleici

Biologie Celulara si Moleculara

Modul Biologie celulara:

Notiuni microscopie

Evaluare organite celulare;

Culturi de celule;

Modul Biologie moleculara:

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Amplificarea ADN/ARN

Electroforeza

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ASO.

RFLP.

DGGE.

Real-Time PCR.


Secvențiere.

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Metode de izolare a acizilor nucleici. Principii.

Consideratii

- tipurile diferite de probe necesita metode diferite de izolare.
- tipuri diferite de probe: tesuturi umane (lichide sau solide), bacterii, plante.
- tipuri diferite de probe: proaspat recoltate, pastrate in diferite conditii (congelate, fixate in parafina), degradate (ex. mumii)
- tipuri diferite de protocoale de lucru: manual, automat
- tipuri diferite de protocoale de lucru in functie de principiul metodei
- sunt similitudini intre protocoale izolare ADN si ARN
- conditii obligatorii pentru a utiliza laboratorul.
- urmarirea unui protocol de lucru

Patologie

- Afectiuni digestive maligne (gastrice, hepatice, pancreatice, colorectale)
- Patologie infectioasa (virus hepatitic B, virus hepatitic C, virusul papiloma uman)
- Ginecologie
- Nefrologie
- etc...

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


The image represents the general outline of how DNA can be extracted from the cell.

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


spalare ADN

rehidratate

DNA extraction step	Chemical
Lysis of cell wall/ cell membrane and Lysis of nuclear membrane	Tris, MgCl ₂ , EDTA, NaCl, SDS, CTAB, Triton X100
Digestion of protein	CTAB, SDS, phenol, chloroform, urea, guanidium isothiocyanate, guanidium thiocyanate, N-Lauroyl sarcosine
Precipitation of DNA	Isopropanol, ethanol, methanol, NaCl, sodium acetate
Washing of DNA	Any alcohol (ethanol commonly used)
Dissolving DNA	TE buffer, distil water

Etapele izolarii ADN

Exemplu de procedeu manual pe baza de solutii

liza celulara


liza nucleara

precipitare proteica

precipitare ADN in izopropanol

spalare ADN cu etanol, uscare

rehidratare cu apa ultrapure sau solutie tampon TE

Etapele izolarii ADN

Exemplu de procedeu manual pe baza de solutii

liza celulara

liza nucleara

precipitare proteica

precipitare ADN in izopropanol

spalare ADN cu etanol, uscare

rehidratare cu apa ultrapure sau solutie tampon TE

Wizard® Genomic DNA Purification Kit

INSTRUCTIONS FOR USE OF PRODUCTS A1120, A1123, A1125 AND A1620.

Quick
PROTOCOL

Isolation of Genomic DNA from Whole Blood

Sample Size	Lysis Solution		Protein Precipitation		DNA Rehydration
	Cell	Nuclei	Solution	Isopropanol	Solution
300µl	900µl	300µl	100µl	300µl	100µl
1ml	3ml	1ml	330µl	1ml	150µl
3ml	9ml	3ml	1ml	3ml	250µl
10ml	30ml	10ml	3.3ml	10ml	800µl

As little as 20µl can be processed with this system. Please see Technical Manual #TM050, Section 3.C.

Red Blood Cell Lysis

- Using volumes from the table above, combine the appropriate volumes of Cell Lysis Solution and blood. Mix by inversion.
- Incubate for 10 minutes at room temperature.
- Centrifuge:
 - ≤300µl sample 13,000–16,000 × g*; 20 seconds
 - 1–10ml sample 2,000 × g; 10 minutes
- Discard supernatant. Vortex pellet.

Nuclei Lysis and Protein Precipitation

- Using volumes from the table above, add Nuclei Lysis Solution and mix by inversion.
- Add Protein Precipitation Solution; vortex for 20 seconds.
- Centrifuge:
 - ≤300µl sample 13,000–16,000 × g*; 3 minutes
 - 1–10ml sample 2,000 × g; 10 minutes

DNA Precipitation and Rehydration


- Transfer supernatant to a new tube containing isopropanol (using volumes from table above). Mix.
- Centrifuge:
 - ≤300µl sample 13,000–16,000 × g*; 1 minute
 - 1–10ml sample 2,000 × g; 1 minute
- Discard supernatant. Add 70% ethanol (same volume as isopropanol).
- Centrifuge as in Step 9.
- Aspirate the ethanol and air-dry the pellet (10–15 minutes).
- Rehydrate the DNA in the appropriate volume of DNA Rehydration Solution for 1 hour at 65°C or overnight at 4°C.

*Maximum speed on a microcentrifuge.

Additional protocol information is available in Technical Manual #TM050, available online at: www.promega.com

ORDERING/TECHNICAL INFORMATION:

www.promega.com • Phone 608-274-4330 or 800-356-9526 • Fax 608-277-2601

Promega

Printed in USA. Revised 4/20
Part #98-B022

Etapele izolarii ADN

Exemplu de procedeu manual pe baza de solutii

liza celulara

liza nucleara

precipitare proteica

precipitare ADN in izopropanol

spalare ADN cu etanol, uscare

rehidratare cu apa ultrapure sau solutie tampon TE

Wizard® Genomic DNA Purification Kit

INSTRUCTIONS FOR USE OF PRODUCTS A1120, A1123, A1125 AND A1620.

Quick
PROTOCOL

Isolation of Genomic DNA from Animal Tissue and Tissue Culture Cells

Prepare Tissues

Tissue Culture Cells: Centrifuge at 13,000–16,000 × g^* for 10 seconds. Wash the cell pellet with PBS, vortex and then add 600µl of Nuclei Lysis Solution and mix by pipetting.

Animal Tissue: Add 10–20mg of fresh or thawed tissue to 600µl of chilled Nuclei Lysis Solution and homogenize for 10 seconds. Alternatively, use 10–20mg of ground tissue. Incubate at 65°C for 15–30 minutes.

Mouse Tail: Add 600µl of chilled EDTA/Nuclei Lysis Solution to 0.5–1cm of fresh or thawed mouse tail. Add 17.5µl of 20mg/ml Proteinase K and incubate overnight at 55°C with gentle shaking.

Lysis and Protein Precipitation

1. Add 3µl of RNase Solution to the cell or animal tissue nuclei lysate and mix. Incubate for 15–30 minutes at 37°C. Cool to room temperature.
2. Add 200µl of Protein Precipitation Solution. Vortex and chill on ice for 5 minutes.
3. Centrifuge at 13,000–16,000 × g^* for 4 minutes.

DNA Precipitation and Rehydration

4. Transfer supernatant to a fresh tube containing 600µl of room temperature isopropanol.
5. Mix gently by inversion.
6. Centrifuge at 13,000–16,000 × g^* for 1 minute.
7. Remove supernatant and add 600µl of room temperature 70% ethanol. Mix.
8. Centrifuge as in Step 6.
9. Aspirate the ethanol and air-dry the pellet for 15 minutes.
10. Rehydrate the DNA in 100µl of DNA Rehydration Solution for 1 hour at 65°C or overnight at 4°C.

*Maximum speed on a microcentrifuge.

Additional protocol information is available in Technical Manual #TM050, available online at: www.promega.com

Cells or tissue
with Nuclei
Lysis Solution.

Add RNase
Solution.
Incubate at 37°C.

Add Protein
Precipitation
Solution.

Centrifuge.

Transfer
supernatant
to new tube
containing
isopropanol.

Centrifuge.

Discard
supernatant.
Add ethanol.

Centrifuge.

Aspirate ethanol.
Air-dry pellet.
Rehydrate DNA.

11000111_1,3A

ORDERING/TECHNICAL INFORMATION:

www.promega.com • Phone 608-274-4330 or 800-356-9526 • Fax 608-277-2601

Promega

Printed in USA, Revised 4/20
Part #918022

Etapele izolarii ADN

Exemplu de procedeu manual pe baza de solutii

liza celulara

liza nucleara

precipitare proteica

precipitare ADN in izopropanol

spalare ADN cu etanol, uscare

rehidratare cu apa ultrapure sau solutie tampon TE

Wizard® Genomic DNA Purification Kit

INSTRUCTIONS FOR USE OF PRODUCTS A1120, A1123, A1125 AND A1620.

Quick
PROTOCOL

Isolation of Genomic DNA from Yeast Cultures or Plant Tissue

Prepare Yeast Lysate

1. Pellet cells from 1ml of culture by centrifugation at $13,000\text{--}16,000 \times g^*$ for 2 minutes.
2. Suspend the cell pellet in 293 μ l of 50mM EDTA.
3. Add 7.5 μ l of 75 units/ μ l lyticase and mix gently.
4. Incubate for 30–60 minutes at 37°C. Cool to room temperature.
5. Centrifuge as in Step 1. Discard the supernatant.
6. Add 300 μ l of Nuclei Lysis Solution. Proceed to **Protein Precipitation and DNA Rehydration**, Step 1 (below).

Prepare Plant Lysate

1. Grind approximately 40mg of leaf tissue in liquid nitrogen.
2. Add 600 μ l of Nuclei Lysis Solution. Incubate at 65°C for 15 minutes.
3. Add 3 μ l of RNase Solution. Incubate at 37°C for 15 minutes. Cool sample to room temperature for 5 minutes. Proceed to **Protein Precipitation and DNA Rehydration**, Step 1 (below).

Protein Precipitation and DNA Rehydration


	Yeast	Plant
1. Add Protein Precipitation Solution. Vortex. For yeast only: Incubate 5 minutes on ice.	100 μ l	200 μ l
2. Centrifuge at $13,000\text{--}16,000 \times g^*$.	3 minutes	3 minutes
3. Transfer supernatant to clean tube containing room temperature isopropanol.	300 μ l	600 μ l
4. Mix by inversion and centrifuge at $13,000\text{--}16,000 \times g^*$.	2 minutes	1 minute
5. Decant supernatant and add room temperature 70% ethanol.	300 μ l	600 μ l
6. Centrifuge at $13,000\text{--}16,000 \times g^*$.	2 minutes	1 minute
7. Aspirate the ethanol and air-dry the pellet.		
8. Add DNA Rehydration Solution.	50 μ l	100 μ l
9. For yeast only: Add RNase. Incubate at 37°C for 15 minutes.	1.5 μ l	—
10. Rehydrate at 65°C for 1 hour or overnight at 4°C.		

*Maximum speed on a microcentrifuge.

Additional protocol information is available in Technical Manual #TM050, available online at: www.promega.com

ORDERING/TECHNICAL INFORMATION:

www.promega.com • Phone 608-274-4330 or 800-356-9526 • Fax 608-277-2601

Promega

Printed in USA. Revised 4/20
Part #9F8022

Etapele izolarii ADN

Exemplu de procedeu manual pe baza de solutii

liza celulara

liza nucleara


precipitare proteica

precipitare ADN in izopropanol

spalare ADN cu etanol, uscare

rehidratare cu apa ultrapure sau
solutie tampon TE

The general outline of DNA extraction to PCR.

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© Genetic Education Inc.